GM501 VitriStore

VitriStore Freeze Kit
Material included in the package:
Product code: 4VF_KIT1

- 1 x 5 ml VitriStore Pre-vitrification Medium (4 VPI-005)
- 1 x 1 ml VitriStore Thaw Medium 1 (4 VT1001)
- 1 x 1 ml VitriStore Freeze Medium 1 (4 VFZ001)
- 1 x 1 ml VitriStore Thaw Medium 4 (4 VT4001)

VitriStore Thaw Kit
Material included in the package:
Product code: 4VT_KIT1

- 1 x 5 ml VitriStore Thaw Medium 1 (4 VT1001)
- 1 x 1 ml VitriStore Thaw Medium 2 (4 VT1002)
- 1 x 1 ml VitriStore Thaw Medium 3 (4 VT1003)
- 1 x 1 ml VitriStore Thaw Medium 4 (4 VT4001)

The media should be used in the order displayed above.

Material not included in the package:

- 4-well dishes
- Freezing tank with liquid nitrogen
- Water bath (able to hold 37°C)
- Pipettes
- Forceps
- Vitrication device
- Water bath (ISO 5 environment)
- Microscope
- Lab timer

Product specifications and quality control:

- All raw materials are of highest available purity (European Pharmacopoeia and/or USP standard) if applicable.
- A certificate of analysis is available for each batch upon request from our website with respective lot number.
- The MSDS for the media are available upon request and can also be downloaded from our website.
- GM501 VitriStore Freeze/VitriStore Thaw are manufactured and tested according to the following specifications:

pH	7.20-7.40
Osmolality (mOsm/kg)	270-290 (VitriStore Freeze PV) 805-850 (VitriStore Thaw 1) 435-565 (VitriStore Thaw 2) 405-435 (VitriStore Thaw 3) 270-290 (VitriStore Thaw 4)
Endotoxins (EU/ml)	< 0.25
Mouse Embryo Assay (MEA) (1-cell assay, blastocysts after 96 h in %)	≥80

EN: Intended use/intended users:

- GM501 VitriStore Freeze/VitriStore Thaw are a set of ready-to-use media for vitrification and thawing of human embryos.
- The intended users are IVF professionals (lab technicians, embryologists or medical doctors).

Components:

- GM501 VitriStore Freeze/VitriStore Thaw are DMSO/Ethylene glycol based vitrification media that also contain PBS, sorbitol, Ficolin and human serum albumin (10-20 g/liter).
- GM501 VitriStore Freeze/VitriStore Thaw do not contain antibiotics.

Instructions for use:

- Instructions are well mixed before use. We strongly advise to read through all the steps of the vitrification/warming procedure before starting the procedure.

Preliminary steps

- In a 4-well dish fill the first well with 300 µl of VitriStore Pre-vitrification Medium, the second with 250 µl VitriStore Freeze Medium 1 and the third with 250 µl VitriStore Freeze Medium 2. Next open as many packs of vitrification devices as will be required for the vitrification step. Conveniently place the separate parts of the vitrification device on the workbench for easy access later in the procedure.

- Use 5 vitrification cycles (of the same patient) can be performed with one media set-up. Do not use the same media for different patients!

Freezing preparation

- Transfer embryos from the blastocyst cell culture medium into each of the VitriStore Freeze solutions using the following scheme:

Stage	VitriStore Pre-vitrification	VitriStore Freeze 1	VitriStore Freeze 2	Temperature
Early Blastocyst, Morulae	2 min.	2 min.	30 sec.	Room temperature
Blastocyst-expanded BI	2 min.	3 min.	30 sec.	37°C
Blastocyst-expanded BI + artificial shrinkage	2 min.	2 min.	30 sec.	Room temperature

**Before starting the vitrification procedure, in order to reduce the negative effect of the blastocyst, expanded blastocysts should be collapsed by reducing artificially with a glass pipette the volume of the blastocyst. (Vanderzwalmen et al., 2002; Son et al., 2003; Hirako, 2004)*

- Using an attenuated pipette or an equally suitable device, deposit maximum 2 blastocysts in a volume of approximately 0.3 µl of VitriStore Freeze 2 in the gutter of the tip of the vitrification straw. Place the vitrification straw in the outer sheath and seal as indicated in the instructions for use of the vitrification device.
- Plunge the sealed device into the liquid nitrogen.

Thawing

- In a 4-well dish, fill the first well with 1 ml VitriStore Thaw Medium 1, the second with VitriStore Thaw Medium 2 (250-300 µl), the third with VitriStore Thaw Medium 3 (250-300 µl), the fourth with VitriStore Thaw Medium 4 (250-300 µl). Preheat the media to 37°C.
- Remove the vitrification straw from the outer sheath as indicated in the instructions for use of the vitrification device. Immediately plunge the vitrification straw into pre-heated VitriStore Thaw Medium 1 (37°C) and leave in Thawing 1 Medium for 3 minutes.
- Transfer into VitriStore Thaw Medium 2 (37°C) and leave in this medium for 2 minutes.
- Transfer into VitriStore Thaw Medium 3 (37°C) and leave in this medium for 2 minutes.
- Transfer into VitriStore Thaw Medium 4 (37°C) and wash for at least 1 minute.
- Transfer into blastocyst culture medium for continued cell culture.

GM501 VitriStore Freeze/VitriStore Thaw and Embryoniculture:

- GM501 VitriStore Freeze/VitriStore Thaw can be used in combination with GM501 Media (Washing- and IVF-Media).

Precautions and warnings:

- Standard measures to prevent infections resulting from the use of medicinal products derived from human donors and plasma products for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of all viruses, including HIV, are included in the instructions for use of the media. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infectious agents cannot totally be excluded. This also applies to unknown or emerging viruses and other pathogens.

- Prevention of HIV or hepatitis: Always wear protective clothing when handling specimens.
- Always work under strict hygienic conditions (LAF-bench ISO class 5) to avoid possible contamination.
- Only for the intended use.
- The long term safety of embryo vitrification on children born following this procedure is unknown.
- The facility of the user is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

Pre-use checks:

- Do not use the product if it becomes discoloured, cloudy, or shows any evidence of microbial contamination.
- Do not use the product if seal of the container is opened or damaged when the product is used.

Storage instructions and stability:

- The shelf life is 12 months from time of manufacture.
- Store between 2-8°C.
- Do not freeze, thaw, or freeze again.
- Keep away from (sun) light.
- The product can be used safely up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8°C.
- Stable after transport (max. 5 days) at elevated temperature (≤ 37°C).
- Content cannot be re-sterilized after opening.
- Do not use after expiry date.

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