# GM501 Hyaluronidase

**Material included in the package:**

- Product code: 4 HY 0001-5
- 5 x 1 ml GM501 Hyaluronidase
- Product code: 4 HY 0101
- 1 x 10 ml GM501 Hyaluronidase

**Material not included in the package:**

- Incubator at 37°C (no CO<sub>2</sub>)
- Plastic dishes
- LAf bench (ISO 5 environment)
- Pipette suitable for denudation of oocytes
- Mineral Oil (e.g. GM501 Mineral Oil)
- Washing medium (e.g. GM501 Wash)
- Microscope

**Product specifications and quality control:**

- All raw materials are of highest available purity (European Pharmacopoeia and/or USP standard) if applicable.
- A certificate of analysis is available for each batch upon request from our website with respect to lot number.
- The MSDS for GM501 Hyaluronidase is available upon request and can also be downloaded from our website.
- GM501 Hyaluronidase is manufactured and tested according to the following specifications:

pH	7.30-7.60
Osmolality (mOsm/kg)	270-290
Sterility	sterile - SAL10% (Sterility Assurance Level)
Endotoxins (EU/ml)	< 1.00
Mouse Embryo Assay (MEA)	≥80 (1-cell assay, blastocysts after 96h in %; 30 min exposure at zygote stage)

**EN: Intended use/intended users:**

GM501 Hyaluronidase is a ready-to-use solution designed to facilitate the mechanical removal of cumulus cells. Hyaluronidase digests the extracellular matrix in the cumulus-oocyte complex consisting of hyaluronic acid.

The intended users are IVF professionals (lab technicians, embryologists or medical doctors).

**Composition:**

- 80 IU/ml pharmaceutical grade hyaluronidase from bovine origin solved in HEPEs-buffered medium, containing 4.00 g/liter human serum albumin.

**Instructions for use:**

GM501 Hyaluronidase contains HEPEs; no CO<sub>2</sub> incubation is required. Just warm it up to 37°C (incubation in a CO<sub>2</sub> incubator will keep the pH below 7).

Depending on the dishes, 100 µl (microdrops) or 400-600 µl (center-IV well) of GM501 Hyaluronidase should be used for denudation. Additionally prepare 3-5 drops (100 µl) of washing medium (e.g. GM501 Wash) for each washing.

The dishes should be covered with Mineral Oil (e.g. GM501 Mineral Oil).

Place the oocyte in the Hyaluronidase (up to max. 5 oocytes) for about 30 seconds.

Using a fine pipette, transfer the partially denuded oocytes to the first washing drop.

- Remove the corona by repeatedly pipetting the oocytes.
- Use the other drops to further wash the denuded oocyte.

**Precautions and warnings:**

Standard measures to prevent infections resulting from the use of potential products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot totally be excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin, manufactured to European Pharmacopoeia specifications by established processes.

Handle all specimens as if capable of transmitting HIV or hepatitis. Always wear protective clothing when handling specimens.

Product does not contain any antibiotics. Always work under strict hygienic conditions (e.g. LAf-bench, ISO Class 5) to avoid possible contamination.

Only for the intended use.

The use of any other device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

**Pre-use checks:**

- Do not use the product if it becomes discoloured, cloudy, or shows any evidence of microbial contamination.
- Do not use the product if the seal of the container is opened or defective when the product is delivered.

**Bovine sourced Hyaluronidase:**

The Pharmaceutical grade hyaluronidase used in this product is sourced from bovine testes. It is certified not to contain any neuronal tissue (Certificate of Suitability for use in humans).

Hyaluronidase is sourced from animals determined „fit for human consumption“ and originating from countries with „negligible BSE risk“ (Resolution No. 18 „Recognition of the Bovine Spongiform Encephalopathy Risk Status of Member Countries (2014)“, adopted by the OIE). According to the OIE guidelines, Tissue Infectivity Distribution in Transmissible Spongiform Encephalopathies (2010), lesions from bovine source are classified as Category IC material (i.e. „Tissues with no detected infectivity“).

**Storage instructions and stability:**

- The shelf life is 12 months from date of manufacture.
- Store at 2-8°C.
- Do not freeze before use.
- Keep away from (sun) light.
- The product can be used safely up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8°C.
- Always after transport (max. 5 days) at elevated temperature (>37°C).
- Do not use after expiry date.

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