



1.5 ml Plain Skirted Tube, Natural Tethered Screw Cap (Sterile)

SKU E1415-2242



Product Highlights

- Routinely tested to 20,000 x g
- Sterile and certified RNase, DNase, DNA, PCR Inhibitor and Endotoxin free.
- Tubes supplied ready assembled with natural tethered caps with o-ring seal.
- O-ring forms a seal against the rim of the tube
- Tube made of polypropylene (PP)

Product Applications

Storage and transportation of biological material

Ideal for boiling and freezing applications

Centrifugation

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

STARLAB International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Pack Size	Bag
Cap style	Natural tethered cap
Tube Style	Plain
Tube Base	Skirted
g-Force max.	20,000 x g
Free of endotoxins	Yes
DNA free	Yes
PCR inhibitor free	Yes
RNase free	Yes
DNase free	Yes
Color	Natural
Autoclavable	No
Sterile	Yes
Volume	1.5 ml
Packaging Item	Tubes
Volume range	1.5 ml
Packaging Unit	Bags
Pack Size	250 Tubes (5 Bags × 50 Tubes)

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

STARLAB International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about 1.5 ml Plain Skirted Tube, Natural Tethered Screw Cap (Sterile)

Assembled sterile tubes & caps

A range of sterile and pre-assembled microtubes, available with a range of different cap styles. Supplied in convenient, resealable mini bags which stand up on your bench. Suitable for many applications from sample storage in the gas phase of liquid nitrogen through freezer to boiling.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

STARLAB International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de