



8-Strip PCR Tubes with Individual Attached Caps, 0.1 ml (natural, Xtra-Clear)

SKU A1402-3800



Product Highlights


- Polypropylene PCR strips
- Attached caps eliminates the possibility of cross-contamination through accidental tube–cap mix up
- Tubes with individually attached caps can easily be split into smaller sections by snapping the strips apart
- Suitable for Realtime PCR (quantitative PCR = qPCR)
- Certified RNase, DNase and DNA free
- Supplied in resealable bags

Product Applications

Standard PCR (Endpoint PCR)

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

STARLAB International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de




General Data

Pack Size	Pack
Cap style	Individually attached flat caps – (Xtra-Clear), low-profile
No. of tubes	8
Material	B-7425 Polypropylene
Free of endotoxins	Yes
DNA free	Yes
PCR inhibitor free	Yes
RNase free	Yes
DNase free	Yes
Color	Natural
Autoclavable	Yes
Sterile	No
Volume	0.1 ml
Packaging Item	Tubes
Volume range	0.1 ml
Packaging Unit	Packs
Pack Size	1 Packs × 120 Tubes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

STARLAB International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about 8-Strip PCR Tubes with Individual Attached Caps, 0.1 ml (natural, Xtra-Clear)

STARLAB PCR strips are manufactured from pure polypropylene. The thin-walled tubes ensure efficient and homogenous heat transfer. Available in 8- or 12-Strip Tubes and Caps.

Low-Profile Strips (0.1 ml): 15.24 mm high, for use in cyclers with low-profile capability only.

Standard Strips (0.2 ml): 20.83 mm high, for use in the majority of cyclers.