

starlab Mixer HC

- mixing: 300 - 1500 rpm
- heating/cooling: 13 - 99 °C
- timing: 00:01 - 99:59 h



9 thermoblocks: MAXIMUM FLEXIBILITY



- 0.5 ml Microcentrifuge Tubes (x24)
- 1.5 ml Microcentrifuge Tubes (x24)
- 2.0 ml Microcentrifuge Tubes (round)(x24)
- 2.0 ml Microcentrifuge Tubes (v-shaped, Ø 12mm)(x24)
- 5.0 ml Preparation Tubes (x8)
- 15 ml Conical Centrifuge Tubes (x8)
- 50 ml Conical Centrifuge Tubes (x4)
- 1.5/2.0 ml Cryo Tubes (x24)
- Micro / Deepwell Plates, with Lid

**bloky i přístroj jsou plně kompatibilní s bloky přístroje i přístrojem Eppendorf Thermomixer comfort*

 **LifeM**
LIFE SCIENCE AND MEDICAL SUPPLIES

Life M s.r.o.
Voděradská 2552/16
251 01 Říčany u Prahy
Česká republika
Tel. +420 323 605 454
IČ: 05081904
DIČ: CZ05081904
<https://www.lifem.cz/>

