


AFI RXN-100



DEDICATED PLASTIC BUCKETS

Model	V	R
Max Capacity	4 x 100ml	
Angle	90°	
RPM Max	4 100	4 100
ACR Max	2 932	2 932



Accessories



Ref.	AFI-RXN100
------	------------

Tubes

Capacity (ml)	1,1-1,4		1,3		5		1,6-5		5-7	
∅ x L (mm)	8x82		10x45		12x75		13x75		13x100	
Tubes per adapter / rotor	8 32		13 52		13 52		7 28		10 40	
Adapter										
Hole ∅ x L (mm)	9,5x55		11x40		13,5x58		13,5x58		13,5x80	
Radius (mm)	120		99		131		132		141	
Max Charge (gr)	200		80g		200		200		200	
Ref.	ISA-PRBT1	ISA-PSBT1	ISA-PRBT1.3	ISA-PSBT1.3	ISA-PSBT5-G12	ISA-PRBT5	ISA-PSBT5	ISA-PRBT6	ISA-PSBT6	

Tubes

Capacity (ml)	2,6-4,9		5-7		2,7-3		4,5-5	
∅ x L (mm)	13x81		13x100		11x82		11x108	
Tubes per adapter / rotor	7 28		7 28		8 32		8 32	
Adapter								
Hole ∅ x L (mm)	13,5x58		13,5x80		12,5x55		12,5x80	
Radius (mm)	131		141		120		137	
Max Charge (gr)	200		200		200		200	
Ref.	ISA-PRBT5	ISA-PSBT5	ISA-PRBT6	ISA-PSBT6	ISA-PRBT5ST	ISA-PSBT5ST	ISA-PRBT6ST	ISA-PSBT6ST

Tubes

Capacity (ml)	4-8,5		7-10		15	50	100
∅ x L (mm)	15x91 108		16x100		17x120	29x115	44x100
Tubes per adapter / rotor	4 16		4 16		3 12	1 4	1 4
Adapter							
Hole ∅ x L (mm)	17,5x 80		17,5x 83		18x80	30x70	46x80
Radius (mm)	142		142		156	147	147
Max Charge (gr)	200		200		100	120	200
Ref.	ISA-PRBT10	ISA-PSBT10	ISA-PRBT10	ISA-PSBT10	ISA-PRBT15C	ISA-PRBT50C	ISA-PRBT100