

DAKEWE

HistoFlo™ Slides



Code	Name	Dimensions	Paint tab color	Description	Packaging
7109030A	Dakewe HistoFlo™ microscope slides	25 × 75 mm, thickness: 1.0 - 1.2 mm	Available in 5 colors: White / Pink / Yellow / Blue / Green	Float glass, FDA approved & CE marked. Ideal for laser slide printers.	50 slides/pack, 50 packs/carton
DKWAS001A	Dakewe HistoFlo™ adhesion microscope slides	25 × 75 mm, thickness: 1.0 - 1.2 mm	Available in 5 colors: White / Pink / Yellow / Blue / Green	Float glass, FDA approved & CE marked. Ideal for laser slide printers. Special coating to improve tissue adherence.	50 slides/pack, 50 packs/carton

