


50 ml Centrifuge Tube, Conical, Loose (Sterile)



Don't get in a spin! Use STARLAB centrifuge tubes!

Product Features

- Clear, pre-steriled polypropylene tubes with blue HDPE screw caps
- White writing surface
- Clearly printed graduations
- Certified RNase, DNase, DNA and Pyrogen free
- Non-Cytotoxic
- Centrifugation Rate Conical: Swing Bucket Rotor: 3,200 x g (swing bucket) / 9,500 x g (fixed angle rotors)

50 ml Centrifuge Tube, Conical, Loose (Sterile)			
	Description	Pack Size	Cat. No.
	50 ml Centrifuge Tube, Conical, Loose (Sterile)	500 Tubes (10 Bags x 50 Tubes)	E1450-0200

Technical Specifications	
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
Free of heavy metals	No
Human DNA free	Yes
Non-cytotoxic	Yes
Pack Size	500 Tubes (10 Bags x 50 Tubes)
Pyrogen free	Yes
RNase free	Yes
Sterile	Yes
Tube Base	Conical
g-Force max.	3,200 x g (swing bucket) / 9,500 x g (fixed angle rotors)
Volume	50 ml