

8-Strip PCR Tubes, Attached Flat Cap Strip, 0.2 ml (Xtra-Clear)




STARLAB PCR strips are manufactured from pure polypropylene. The thin-walled tubes ensure efficient and homogenous heat transfer. Available in 8- or 12-Strip Tubes and Caps.

Product Features

- Polypropylene PCR strips
- Attached caps eliminates the possibility of cross-contamination through accidental tube-cap mix up
- Tubes with individually attached caps can easily be split into smaller sections by snapping the strips apart
- Suitable for Realtime PCR (quantitative PCR = qPCR)
- Certified RNase, DNase, DNA and PCR Inhibitor free
- Autoclavable
- Supplied in resealable bags

8-Strip PCR Tubes, Attached Flat Cap Strip, 0.2 ml (Xtra-Clear)

	Description	Pack Size	Cat. No.
	8-Strip PCR Tubes, Attached Flat Cap Strip, 0.2 ml (Xtra-Clear)	1 Packs × 125 Tubes	I1402-2700

Technical Specifications

Autoclavable	Yes
Cap style	Attached flat cap strip (Xtra-Clear)
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
Material	B-7425 Polypropylene
No. of tubes	8
PCR inhibitor free	Yes
Pack Size	1 Packs × 125 Tubes
RNase free	Yes
Sterile	No
Color	Natural
Volume	0.2 ml