
20 μ l UltraPoint[®] TipOne Graduated Filter Tip (sterile), Racked



Nothing goes through!

Protect your samples against unwanted contamination by using TipOne® filter tips. The TipOne® filter is one of the best on the market and protects against aerosol cross-contamination without trapping your valuable samples. The filters contain no additives to interfere with samples.

During pipetting it is not possible to see with the eye if fine aerosols are going into your pipette. If they are, those aerosols can contaminate other samples at a later date.


TipOne® filter tips are proven to block >99% of aerosols and are among the best filters on the market. The results, shown right, speak for themselves!

**The comparative tests were conducted in an independent test institute with filter tips from other leading manufacturers.*

Product Features

- **Graduated** at 10 µl, 50 µl and 100 µl
- **Quality:** Highly effective HDPE filter
- **Safe:** Filter contains no additives to interfere with your samples
- **Retrieval:** The filter will not absorb samples, allowing recovery if necessary
- **Certified:** Sterile, and RNase, DNase, DNA and Pyrogen free
- **Sterile**
- **UltraPoint®** thin-wall orifice helps prevent sample hang up and makes it easier to touch off the last drop

20 µl UltraPoint® TipOne Graduated Filter Tip (sterile), Racked

	Description	Pack Size	Cat. No.
	20 µl UltraPoint® TipOne Graduated Filter Tip (sterile), Racked	960 Tips (10 Racks × 96 Tips)	S1123-1810

Technical Specifications

Autoclavable	No
Color Code Wafer	Yellow
Material	Polypropylene
Pack Size	960 Tips (10 Racks × 96 Tips)
Sterile	Yes
Color	Natural
Volume	20 µl
For use with	Abimed®, Biohit®, BRAND®, Capp®, Eppendorf®, Gilson®, HTL, Nichiryo®, Rainin®, Sartorius®, Socorex®, STARLAB®, Thermo Scientific™