

Racking guide: upright freezers.

| | EPENDORF (NEW BRUNSWICK) | | | | | | | PHCBI (PANASONIC / SANYO) - MDF RANGE | | | | | | | | | | HAIER | | | | | | THERMO |
|--|--------------------------|--|--|--|--|--|--|---------------------------------------|--|--|--|--|--|--|--|--|--|-------|--|--|--|--|--|--------|
|--|--------------------------|--|--|--|--|--|--|---------------------------------------|--|--|--|--|--|--|--|--|--|-------|--|--|--|--|--|--------|




| | | | F570h CryoCube | F740h CryoCube | U360VIP | U535VIP | U570 | U725VIP | DU302VX | DU502VH-PH DU502VX | DU702VH-PE DU702VX | DU900V-PE | U33V | U4186S | U500VX | U53V | U5386S | U55V | U700VX | U7386S | U76V | DW86L288 | DW86L338(J) | DW86L388(A/J) | DW86L486 | DW86L490(J) | DW86L628 | DW86L728(J) | ExF60086 | | | | | | |
|---|--------------------|--------------------|---|-------------------|---------|---------|------|---------|---------|-----------------------|-----------------------|-----------|------|--------|--------|------|--------|------|--------|--------|------|----------|-------------|---------------|----------|-------------|----------|-------------|----------|----|--|--|--|--|--|
| SIDE ACCESS FOR STANDARD STORAGE BOXES up to 136 x 136 x 53 mm | | |  | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| E6000-0033 | 9-place (3d x 3h) | 139 x 426 x 169 mm | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| E6000-0034 | 12-place (3d x 4h) | 139 x 426 x 224 mm | | | | | | | | | | | | 8* | | | | | | | | | | | | | | | | | | | | | |
| E6000-0035 | 15-place (3d x 5h) | 139 x 426 x 297 mm | | | | | | | | | | | | 8* | | | | | | | | | 12 | | 16 | | | | | | | | | | |
| E6000-0043 | 12-place (4d x 3h) | 139 x 565 x 169 mm | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| E6000-0044 | 16-place (4d x 4h) | 139 x 565 x 224 mm | 25 | | | | 25 | | | | | | 6* | | | | 16 | | | | | | | | | | | | | | | | | | |
| E6000-0045 | 20-place (4d x 5h) | 139 x 565 x 279 mm | | | | | | | 12 | | | 14* | 6* | | 8* | 8* | | 8* | 12* | 24 | 12* | | 6* | | | 8* | 20 | 24 | | | | | | | |
| E6000-0046 | 24-place (4d x 6h) | 139 x 565 x 334 mm | | | | | | | | 16 | 24 | 14* | | | 8* | 8* | | 8* | 12* | | 12* | | | | | | | | | | | | | | |
| PULL OUT SHELVES FOR STANDARD STORAGE BOXES up to 136 x 136 x 53 mm | | |  | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| E6001-0033 | 9-place (3d x 3h) | 140 x 450 x 173 mm | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| E6001-0034 | 12-place (3d x 4h) | 140 x 450 x 231 mm | | | | | | | | | | | | 8* | | | | | | | | | | | | | | | | | | | | | |
| E6001-0035 | 15-place (3d x 5h) | 140 x 450 x 288 mm | | | | | | | | | | | | 8* | | | | | | | | | 12 | | 16 | | | | | | | | | | |
| E6001-0043 | 12-place (4d x 3h) | 140 x 569 x 174 mm | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| E6001-0044 | 16-place (4d x 4h) | 140 x 569 x 231 mm | 25 | | | | 25 | | | | | | 6* | | | | 16 | | | | | | | | | | | | | | | | | | |
| E6001-0045 | 20-place (4d x 5h) | 140 x 569 x 281 mm | | | | | | | 12 | | | 14* | 6* | | 8* | 8* | | 8* | 12* | 24 | 12* | | 6* | | | 8* | 20 | 24 | | | | | | | |
| E6001-0046 | 24-place (4d x 6h) | 140 x 569 x 337 mm | | | | | | | | 16 | 24 | 14* | | | 8* | 8* | | 8* | 12* | | 12* | | | | | | | | | | | | | | |
| E6001-0047 | 28-place (4d x 7h) | 140 x 569 x 403 mm | | 18 | 9 | 12 | | 18 | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| PULL OUT SHELVES FOR STARSTORE 100 STORAGE BOXES 140 x 155 x 55 mm | | |  | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| E6001-0133 | 9-place (3d x 3h) | 147 x 500 x 179 mm | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| E6001-0134 | 12-place (3d x 4h) | 147 x 500 x 239 mm | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| E6001-0144 | 16-place (4d x 4h) | 147 x 655 x 239 mm | | | | | | | | | | | | | | | | | | | | | | | | | | | | 30 | | | | | |

Table shows the maximum number of racks that fit into the model of freezer.

*These freezers hold racks of different sizes on different shelves; table shows maximum number of each rack size. The maximum number of racks per freezer is the total of each added together.

Version 1: more info to be added as soon as it becomes available.

 starlab.click/uprighttracks

Racking guide: chest freezers.



RACKING FOR STANDARD STORAGE BOXES
up to 136 x 136 x 53 mm

| | | | EPPENDORF (NEW BRUNSWICK) | | | | PHCBI (PANASONIC / SANYO) - MDF RANGE | | | | ARTICO | THERMO |
|------------|---------|--------------------|---------------------------|---------|------|---------|---------------------------------------|-----|-----|-----|----------|-------------|
| | | | C340 | C585VIP | C660 | C760VIP | C2156-VAN | 394 | 594 | 794 | ULTF 420 | CXF ULT2090 |
| E6010-0009 | 9-high | 141 x 148 x 501 mm | | | | | | 21 | | | | |
| E6010-0010 | 10-high | 141 x 148 x 556 mm | | | | | 14 | | | | | |
| E6010-0011 | 11-high | 141 x 148 x 611 mm | | | | | | | | 24 | | |
| E6010-0012 | 12-high | 141 x 148 x 666 mm | | | | | | | | | | 33 |
| E6010-0013 | 13-high | 141 x 148 x 721 mm | 13 | 13 | 13 | 13 | | 24 | 36 | | | |

Table shows the maximum number of racks that fit into the model of freezer.

 starlab.click/chestracks

SUSTAINABILITY JUST GOT COOLER!

- Freezer racking helps optimise capacity by keeping freezers organised and maximising storage space for efficient energy use
- Made from 85% recycled stainless steel
- Manufactured in a factory powered by 100% sustainable energy
- Made in Germany. Short journey from factory to customer reduces carbon footprint

#EcolutionMovement

