

# Screw Cap Microtubes.













## THE TUBE RANGE

Starlab's Screw Cap Microtubes consist of a range of tubes and caps which are suitable for many applications from sample storage in the **gas phase of liquid nitrogen** to freezer, through boiling to autoclave. Non-sterile tubes and caps are supplied separately. Sterile tubes and caps are supplied pre-assembled.

**All non-sterile tubes/caps certified RNase, DNase, DNA free. Tested Endotoxin-free.**



### TUBES ONLY

										
<b>Tube Volume:</b>	0.5 ml	0.5 ml	0.5 ml	1.5 ml	1.5 ml	1.5 ml	1.5 ml	2.0 ml	2.0 ml	2.0 ml
<b>Tube Style:</b>	Ribbed	Plain	Plain	Ribbed	Ribbed	Graduated	Plain	Ribbed	Graduated	Graduated
<b>Tube Base:</b>	Skirted	Skirted	Conical	Skirted	Conical	Conical	Skirted	Skirted	Conical	Skirted
 <b>Natural</b>	<b>E1405-2120</b> Pack 500	<b>E1405-2140</b> Pack 500	<b>E1405-2130</b> Pack 500	<b>E1415-2220</b> Pack 500	<b>E1415-2210</b> Pack 500	<b>E1415-2230</b> Pack 500	<b>E1415-2240</b> Pack 500	<b>E1420-2320</b> Pack 500	<b>E1420-2330</b> Pack 500	<b>E1420-2340</b> Pack 500
 <b>Amber</b> for light sensitive samples	<b>E1405-2127</b> Pack 500	<b>E1405-2147</b> Pack 500	n/a	<b>E1415-2227</b> Pack 500	n/a	<b>E1415-2237</b> Pack 500	<b>E1415-2247</b> Pack 500	n/a	<b>E1420-2337</b> Pack 500	<b>E1420-2347</b> Pack 500

### TUBES

- 100% polypropylene
- Routinely tested to 20,000 xg
- Autoclavable
- Amber options for light sensitive samples

### CAPS

- Easy 3/4 turn close/open
- Autoclavable
- Three cap styles available




See page 2 for cap ordering information.

 [starlab.click/screwcaptubes](http://starlab.click/screwcaptubes)

# Screw Cap Microtubes.









## CAP OPTIONS

### EASYGRIP CAPS



Opaque white<sup>1)</sup> polypropylene cap with coloured seal/grip (inert TPE). One-piece construction provides a seal that will never fall out!

<sup>1)</sup>E1480-0300 is not opaque

								
Natural/White	White/Blue	White/Green	White/Orange	White/Red	White/Violet	White/Yellow	White/White	Mixed
<b>E1480-0300</b> Pack 500	<b>E1480-0301</b> Pack 500	<b>E1480-0302</b> Pack 500	<b>E1480-0303</b> Pack 500	<b>E1480-0304</b> Pack 500	<b>E1480-0305</b> Pack 500	<b>E1480-0306</b> Pack 500	<b>E1480-0309</b> Pack 500	<b>E1480-0399</b> Pack 500

 [starlab.click/screwcaps](http://starlab.click/screwcaps)

### CAPS WITH O-RING



**Standard Cap**  
Polypropylene cap with separate, ethylene-propylene o-ring to form a seal against rim of tube.




**Tethered Cap**  
Polypropylene cap with separate, ethylene-propylene o-ring to form a seal against rim of tube.

										
Natural	Blue	Green	Orange	Red	Violet	Yellow	Amber	Black	White	Mixed
<b>E1480-0100</b> Pack 500	<b>E1480-0101</b> Pack 500	<b>E1480-0102</b> Pack 500	<b>E1480-0103</b> Pack 500	<b>E1480-0104</b> Pack 500	<b>E1480-0105</b> Pack 500	<b>E1480-0106</b> Pack 500	<b>E1480-0107</b> Pack 500	<b>E1480-0108</b> Pack 500	<b>E1480-0109</b> Pack 500	<b>E1480-0199</b> Pack 500
<b>E1480-0200</b> Pack 500	<b>E1480-0201</b> Pack 500	<b>E1480-0202</b> Pack 500	<b>E1480-0203</b> Pack 500	<b>E1480-0204</b> Pack 500	<b>E1480-0205</b> Pack 500	<b>E1480-0206</b> Pack 500	<b>E1480-0207</b> Pack 500	<b>E1480-0208</b> Pack 500	<b>E1480-0209</b> Pack 500	<b>E1480-0299</b> Pack 500

# Screw Cap Microtubes.

## STERILE TUBE & CAP ASSEMBLIES

### STERILE TUBES & CAPS




Tube Volume:	0.5 ml	1.5 ml	1.5 ml	2.0 ml	2.0 ml
Tube Style:	Plain	Graduated	Plain	Graduated	Graduated
Tube Base:	Skirted	Conical	Skirted	Conical	Skirted



with Natural Standard Caps

<b>E1405-2141</b> Pack 5 x 50	<b>E1415-2231</b> Pack 5 x 50	<b>E1415-2241</b> Pack 5 x 50	<b>E1420-2331</b> Pack 5 x 50	<b>E1420-2341</b> Pack 5 x 50
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------



with Natural Tethered Caps

<b>E1405-2142</b> Pack 5 x 50	<b>E1415-2232</b> Pack 5 x 50	<b>E1415-2242</b> Pack 5 x 50	<b>E1420-2332</b> Pack 5 x 50	<b>E1420-2342</b> Pack 5 x 50
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------



with Natural/White EasyGrip Caps

<b>E1405-2143</b> Pack 5 x 50	<b>E1415-2233</b> Pack 5 x 50	<b>E1415-2243</b> Pack 5 x 50	n/a	<b>E1420-2343</b> Pack 5 x 50
----------------------------------	----------------------------------	----------------------------------	-----	----------------------------------



with White/Colour EasyGrip Caps<sup>1</sup>

<b>E1405-030#</b> Pack 5 x 50	<b>E1415-030#</b> Pack 5 x 50	<b>E1415-031#</b> Pack 5 x 50	<b>E1420-030#</b> Pack 5 x 50	<b>E1420-031#</b> Pack 5 x 50
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------

<sup>1</sup>To order the required colour of grip/seal required, replace # as last digit of Article No. with:



1 for Blue



2 for Green



4 for Red



6 for Yellow



### PRE-ASSEMBLED STERILE TUBES AND CAPS

- Supplied in convenient resealable mini pouches which stand up on your desk
- Routinely tested to 20,000 xg
- Wide range of tube/cap options

Sterile tube/cap assemblies are certified RNase, DNase, DNA and Endotoxin free.

See page 2 for more details on cap styles.



Sterile product is supplied in convenient, resealable mini pouches



[starlab.click/sterilemicrotubes](http://starlab.click/sterilemicrotubes)