

STRIP-ID

PIPETTE FOR OOCYTES MANIPULATION

The Strip-ID is a single-use IVF laboratory micropipette used for oocytes manipulation and handling of embryos and oocytes

FEATURES

- superior quality that ensures smooth edges for precise denudation
- strict individual microscopic quality inspection
- flexible, unbreakable material that will not scratch dishes
- convenient colour-coded multi-packaging that allows easy release of a single pipette
- extensive variety of sizes available to suit your application
- fits the "Pipet-ID(s)" and is compatible with most popular used micropipettors

ORIGINAL PICTURE



ORIG. PIC. 200X
ENLARGED

GYNO INFO STRIP-ID

Length	90 mm		
Inner diameter μm	125 μm	145 μm	200 μm
	130 μm	150 μm	275 μm
	135 μm	170 μm	300 μm
	140 μm	175 μm	
Sterile	sterilized by gamma irradiation		
Packaging	10 pipettes / vial 15 vials / inner box 75 vials of 10 pipettes / shipping box		
Quality control	MEA and LAL		
Order number	#3000 Strip-ID125/10 #3000 Strip-ID130/10 #3000 Strip-ID135/10 #3000 Strip-ID140/10 #3000 Strip-ID145/10 #3000 Strip-ID150/10 #3000 Strip-ID170/10 #3000 Strip-ID175/10 #3000 Strip-ID200/10 #3000 Strip-ID275/10 #3000 Strip-ID300/10		

GMP Belgium Strip-ID 12/2016



Gynetics Medical Products n.v.
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium
tel. +32 11 64 58 72 / fax +32 11 64 06 77
info@gynetics.be / www.gynetics.be

GYNETICS[®]
Practise better Medicine

GYNO INFO

STRIP-ID

Order number	ID	Colour
#3000 Strip-ID125/10	125 μm	blue
#3000 Strip-ID130/10	130 μm	yellow
#3000 Strip-ID135/10	135 μm	yellow
#3000 Strip-ID140/10	140 μm	white
#3000 Strip-ID145/10	145 μm	white
#3000 Strip-ID150/10	150 μm	black
#3000 Strip-ID170/10	170 μm	red
#3000 Strip-ID175/10	175 μm	red
#3000 Strip-ID200/10	200 μm	dark red
#3000 Strip-ID275/10	275 μm	green
#3000 Strip-ID300/10	300 μm	transparant



Easy release of a single pipette



Colour-coded multi-packaging

