

# Oosafe®

## CLASS IIA MEDICAL DEVICE DISINFECTANTS

Manufactured in accordance  
EN ISO 13485:2016


CE 0068



# Oosafe® CLASS IIA MEDICAL DEVICE DISINFECTANTS

## Oosafe® DISINFECTANT

\*Effective against enveloped viruses including *HIV*, *Hepatitis B* and *C*, *Herpes simplex*, *Influenza A*, *Ebola*.



Product code: **OOMD-00050**  
Oosafe® Disinfectant,  
50 ml in spray bottle



Product code: **OOMD-01000**  
Oosafe® Disinfectant,  
1 Liter bottle with spray



Product code: **OOMD-05000**  
Oosafe® Disinfectant,  
5 Liters refill container



Product code: **OOMDW-120**  
Oosafe® Disinfectant Wipes,  
70 wipes per bucket

### Features

- Ready-to-use solution
- Valid for 3 years, even after opening
- Transparent and odourless

### Health and safety

- MEA and HSSA tested
- Alcohol-free
- Non-VOC releasing

### Intended use:

For disinfection of medical devices - dental and surgical instruments.

### Application:

Ensure that surface is cleaned before applying disinfectant.

**Spray:** Apply undiluted. Spray and distribute evenly on the surface with a sterile cloth. Wait 15 minutes for the disinfection effect.

**Wipes:** Wipe the surface evenly and wait 15 minutes for the disinfection effect.

### ORDER INFORMATION

| Product code |                           |
|--------------|---------------------------|
| OOMD-00050   | 50 ml in spray bottle     |
| OOMD-01000   | 1 Liter bottle with spray |
| OOMD-05000   | 5 Liters refill container |
| OOMDW-120    | 70 wipes per bucket       |

